

IMPORTANT INFORMATION REGARDING THE PATIENT COMPLIANCE STATISTICS BEING RECORDED BY THE RESPLUS

The TF/SD card for the ResPlus Auto CPAP machine comes pre-loaded with the V Health Respiratory Management System App. Quick start guides on how to use the TF card and install the VHealth App are included with each ResPlus unit.

SUMMARY OF RECORDED SLEEP STATISTICS:

USER PROFILE

- Personal information

DATA RECORDING

- Start Date
- End Date

USAGE:

- Total Usage
- Total days
- Total hours used
- Average daily usage
- Used days greater than 4 hours

RESPIRATORY PARAMETERS:

- Apnea index
- Hypopnea Index
- AHI
- ODI

USAGE INFORMATION:

- Pressure – max – average
- Leak – litres/minute – max – average
- Tidal Volume – max – average
- BPM – max – average

DEVICE SETTINGS:

- Working mode
- Humidifier setting
- Start pressure
- Max pressure
- Min pressure
- Ramp set point
- Auto on/Auto off set

Graphs are provided showing - daily usage – AHI – Pressure – leak.
Trend charts are also provided showing daily – AHI – Pressure – leak.
The ability to PRINT any of the above reports/parameters is available.

FOR FULL DETAILS GO TO RESPLUS – V HEALTH – USER MANUAL PDF IN OUR “RESOURCES” SECTION ON OUR WEBSITE.